Ortho-Strips Coarse

DE

Ortho-Strips Coarse

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Ortho-Strips Coarse

University of Zurich

Oscillating instruments for a fast and efficient reduction of interproximal enamel in orthodontic treatments.

To be used in combination with the contra-angles:

- Intensiv Swingle, WG 69 LT Ortho, with light or
- Intensiv Swingle, WG 69 A Ortho, without light

Alternatively:

- KaVo 61 LG
- KaVo contra-angle Intracompact E61C

Indications

Unilateral and bilateral interproximal enamel reduction during orthodontic IPR protocol.

Orthodontic treatments often require the reduction of interproximal teeth dimensions. An adequate interproximal stripping methodology ("InterProximal Reduction" = IPR) involves a controlled and efficient enamel reduction.

The reduction of interproximal enamel is needed to create necessary space to adjust teeth in a correct and harmonious position. After the reduction of enamel by means of a coarse diamond grit-size (80 µm), it is necessary to proceed with the finishing and polishing of surface roughness.

Recommendations for use

- Contra-angle turning approx. 4000 rpm.
- In case of very tight contact points we suggest to separate the teeth by placing a wooden wedge cervically.
- Strip introduction: whilst the contra-angle is running, insert the strip in an oral to buccal cavity direction. This prevents the instruments from jamming.
- Contact pressure: bend the strip 1-1.5 mm to achieve the ideal 1-2.5 Newton contact pressure.
- Avoid bending and levering to a strong extent: use of excessive force (strong levering effect with the contra-angle) could lead to the breakage of the instrument.
- Do not stop the instrument in the interdental space.

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