DE

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EN

# EN. Processing (cleaning, disinfection, and sterilization) of diamond-coated instruments

## 1 Fundamental remarks

All instruments are to be cleaned, disinfected, and sterilized prior to each application; this is required as well prior to the first use after receipt of the non-sterile instruments (cleaning and disinfection after removal of the protective packaging, sterilization after packaging in dedicated sterilization pouches). A sufficient cleaning and disinfection (to be confirmed as described in chapter 7) is an indispensable requirement for an effective sterilization of the instruments.

The dedicated person for processing is in charge for the sterility of the instruments. Therefore, the dedicated person has to ensure that only sufficiently device and product specifically validated procedures will be used for cleaning, disinfection, and sterilization, that the used devices (Washer-Disinfector = WD, sterilizer) will be maintained and checked regularly, as well as that the validated parameters will be applied for each cycle.

Attention is needed to avoid an additional unwanted contamination of the sterilization tray by not placing used instruments back into the tray (separate collection). Pre-cleaning, cleaning, disinfection, and checking of the contaminated instruments is needed. Then the instruments need to be placed back into the sterilization tray and the completely equipped sterilization tray needs to be sterilized.

Additionally, attention is required to the legal provisions valid as well as to the processing procedures of the doctor's practice (clinic or hospital). This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA), which can require the application of cleaning detergents with proven prion efficiency as well as a sterilization with more intensive parameters.

Isolate contaminated material with suitable packaging and labelling.

General precautions should be observed by all medical facility personnel working with contaminated or potentially contaminated medical devices. Care should be taken when handling devices with sharp tips or cutting edges.

Personal protective equipment (PPE) must be worn when handling or working with contaminated or potentially contaminated materials, devices and equipment. The personal protective equipment must comply with the regulations of the respective country and must include at least the following:

- Protective clothing (liquid-tight),
- masks,
- eye and face protection,
- gloves and
- overshoes

Based on the intended use and the geometry/surface properties, we recommend classification as semi-critical B or critical B (depending on the type and location of actual use).

Do not use any fixing cleaning agents or disinfectants.

If unsuitable cleaning agents or disinfectants are used or if unsuitable procedures are applied, this can lead to the loss of warranty claims and have negative consequences for the instruments:

- Damage
- corrosion
- Discoloration of the product
- the necessary hygiene is not guaranteed

## 2 Restriction of reprocessing

### 2.1 Material aspects

Attention is necessary that the listed substances are not ingredients of the cleaning or disinfection detergent:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 12, neutral/ enzymatic or weak alkaline cleaner recommended)
- organic solvents (for example: acetone, ether, alcohol, benzene)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

When selecting the detergents, it has in addition to be considered that corrosion inhibitors, neutralizing agents, and/ or rinse aids may cause potential critical remnants on the instruments.

Acid neutralizing agents or rinse aids must not be applied.

Instruments and sterilization trays must not be cleaned by use of metal brushes or steel wool to avoid scratches.

Instruments and sterilization trays must not be exposed to temperatures higher than 142 °C (288 °F) (see chapter 10).

In case of application of sterilization trays of other manufacturers, attention to the instructions for use and processing of the specific manufacturer is required.

### 2.2 Processing cycles

The number of admitted processing cycles depends on the state of the instruments. High attention needs to be given on adequate care and checking for wearing out state, damage and cleanliness of the instruments (see chapter 7). Used up or damaged instruments should not be further reprocessed or used. The end of the service life (time at which reprocessing can no longer be considered safe) is defined by the use of the instruments. Defective instruments, e.g. with missing diamond coating, fractured working parts, corroded surfaces and/or with detached color markings as well as bent instruments, must be disposed of by the operator and must not be used.

The user is responsible for each further processing as indicated before as well as for the use of damaged and unclean instruments.

Any liability is excluded in the event of non-compliance.

### 3 Preparations at point of use

Not applicable (see pre-cleaning, chapter 5.2)

### 4 Storage and transport to the reprocessing area

To avoid damage, the instruments should be stored safely and transported in a closed container to prevent contamination of the environment until the point of further reprocessing.

Transport the instruments to the reprocessing room and place them in a basin with cleaning solution.

## 5 Cleaning and disinfection

### 5.1 Basics

The pre-cleaning step is to be performed both before a manual or automated cleaning and disinfection procedure.

If available, an automated procedure (Washer-Disinfector = WD) should be used for cleaning and disinfection of the instruments. In case a WD is not available, a manual procedure, inclusive of application of an ultrasonic bath should be used. The significantly lower efficiency and reproducibility of a manual procedure has to be considered.

### 5.2 Pre-cleaning

Equipment:

- ultrasonic bath (35 kHz)
- suitable cleaner
- water (max. 10 germs/ml) and low endotoxin (max. 0.25 endotoxin units/ml)
- soft brush (e.g. nylon brush/toothbrush)

Visible impurities must be removed from the instruments directly after application (within a maximum of 2 h).

Procedure:

1. Rinsing the instruments for at least 1 min under running tap water (temperature < 35 °C/95 °F).
2. Soaking the instruments for the given soaking time in the pre-cleaning solution<sup>1</sup> (ultrasonic bath, ultrasound not activated) so that the instruments are sufficiently covered. Attention is necessary that there is no overlapping of the instruments. Cleaning is to be assisted by careful brushing of all surfaces with a soft brush (at beginning of soaking).
3. Activating ultrasound (35 kHz) for an additional soaking time (but not less than 5 min).
4. Removal of the instruments from the pre-cleaning solution and post-rinsing at least three times intensively (at least 1 min) with water.
5. Checking the instruments for visible remnants. If visible remnants are still present, steps 2 to 4 need to be repeated and another visual check needs to be performed. If remnants are still present after pre-cleaning, the instrument must be disposed of.

<sup>1</sup> In case of application of a cleaning and disinfection detergent for this (e.g. in consequence of personnel's safety), it should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example

NOTES:

Pre-cleaning of the instruments should be carried out immediately after treatment, taking personal protection into account. The aim is to prevent organic material and chemical residues from drying on the instruments and to avoid contamination of the environment.

Attention to following points must be given during selection of the cleaning detergent<sup>1</sup>:

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- compatibility of the cleaning detergent with the instruments (chapter 2)

Attention is necessary to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing. Freshly prepared solutions need to be used, as well as sterile or low contaminated water (max. 10 germs/ml) and low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

A disinfectant used in the pre-cleaning step serves only the personnel's safety but cannot replace the disinfection step later to be performed after cleaning.

### 5.3 Automated cleaning/disinfection (Washer-Disinfector = WD)

Equipment:

- Washer/disinfector
- suitable cleaner
- water (max. 10 germs/ml) and low endotoxin (max. 0.25 endotoxin units/ml)
- small parts basket

Attention to following points during selection of the WD is necessary:

- fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- possibility for an approved program for thermal disinfection (A0 value · 3000 or – in case of older devices – at least 5 min at 90 °C/194 °F;
- in case of chemical disinfection, the danger of remnants of the disinfectant on the instruments has to be considered
- fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- post-rinsing with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- drying with filtered air (oil-free, low contamination with microorganisms and particles)

VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the instruments (chapter 2)

- regular maintenance and checking/calibration of the WD

Attention to following points during selection of the cleaning detergent is necessary:

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- additional application – in case of non-application of a thermal disinfection – of a suitable disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/ clearance/ registration or CE marking) compatible with the used cleaning detergent
- compatibility of the used detergents with the instruments (chapter 2)

Attention is necessary to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

After pre-cleaning acc. point 5.2, transfer the instruments to the washer/disinfector

Procedure:

1. Transferring the instruments by the use of a small parts basket in the WD (attention that the instruments do not overlap in the basket).
2. Starting the program.
3. Removal of the instruments from the WD after end of the program.
4. Checking and packing the instruments for sterilization immediately after the removal (chapters 7, 8 and 9, if necessary after additional post-drying at a clean place).

*The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory. The cleanability was validated by application of an ultrasonic bath Sonorex Super RK 514 H (BANDELIN electronic, Berlin, Germany) using neodisher® Mediclean Forte (0.5%, Fa. Dr. Weigert) for pre-cleaning. A washer-disinfector (WD) type Miele Professional PG8581 (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh Germany) was used for the reprocessing using the cleaning agent neodisher® Mediclean Forte (0.5%, Fa. Dr. Weigert) in the Vario TD program.*

*Step 1 Pre-cleaning with cold tap water for 1 minute.*

*Step 2 Cleaning with 0.5% cleaner 55°C for 5 minutes with demineralised water.*

*Step 3 Rinse with demineralised water for 1 minute.*

*Step 4 Thermal disinfection with demineralised water >90°C for 5 minutes.*

*After thermo disinfection, a drying step can be included in the program of the washer-disinfector.*

*The selected parameters corresponded to the Vario-TD program from Miele.*

## 5.4 Manual cleaning and disinfection

Equipment:

- ultrasonic bath (35 kHz)
- suitable cleaner
- water (max. 10 germs/ml) and low endotoxin (max. 0.25 endotoxin units/ml)
- soft brush (e.g. nylon brush/toothbrush)

Attention to following points during selection of the cleaning and disinfection detergents is necessary:

- fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- application of a disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/ registration or CE marking) compatible with the used cleaning detergent
- compatibility of the used detergents with the instruments (chapter 2)

Combined cleaning/disinfection detergents should not be used.

Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used.

Attention is necessary to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Freshly prepared solutions need to be used, as well as sterile or low contaminated water (max. 10 germs/ml) and low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

Procedure:

After pre-cleaning acc. point 5.2, transfer the instruments to the ultrasonic bath.

Cleaning

1. Soaking the instruments for the given soaking time in the cleaning solution<sup>2</sup> (ultrasonic bath, ultrasound not activated) so that the instruments are sufficiently covered. Attention is necessary that there is no overlapping of the instruments. Cleaning is to be assisted by careful brushing of all surfaces with a soft brush.
2. Activating ultrasound (35 kHz) for an additional soaking time (but not less than 5 min).
3. Removal of the instruments from the cleaning solution and post-rinsing at least three times intensively (at least 1 min) with water.
4. Checking the instruments for remnants

Disinfection

5. Soaking the cleaned and checked instruments for the given soaking time in the disinfectant solution<sup>3</sup> so that the instruments

<sup>2</sup> Cidezyme® (Johnson & Johnson GmbH, Norderstedt) was used as a cleanser at a concentration of 0.8% and a holding time of 5 minutes at 40°C and 35 kHz in a Bandelin RK 514 H ultrasonication device.

<sup>3</sup> Cidex OPA (Johnson & Johnson GmbH, Norderstedt) was used as a disinfectant (undiluted) with a holding time of 5 minutes in a disinfection bath at room temperature (23+3°C).

are sufficiently covered. Attention is necessary that there is no overlapping of the instruments.

6. Removal of the instruments from the disinfectant solution and post-rinsing at least five times intensively (at least 1 min) with water.
7. Drying and packing immediately the instruments for sterilization after the removal (chapters 6 and 9 if necessary after additional post-drying at a clean place).

*The fundamental suitability of the instruments for an effective manual cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§15 (5) MPG) test laboratory by application of the cleaning detergent Cidezyme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson GmbH, Norderstedt) considering the specified procedure.*

**NOTE:**

Chemical disinfection harbors the risk of disinfectant residues on the instruments and can affect the biocompatibility of the products, therefore thorough rinsing after the disinfection step is absolutely necessary.

**6 Drying**

In case that moisture remains on the product after the automated reprocessing and in case that you use a manual process, further manual drying of the products is necessary.

To avoid any kind of deterioration of the instrument, make sure that the instrument is completely dry after each cycle. If moisture remains on or in the instrument after the cleaning process, dry it with a lint-free cloth and drying with medical compressed air if necessary.

**7 Checking**

All instruments need to be visually checked after pre- cleaning and manual or automated cleaning and disinfection, respectively, on corrosion, damaged surfaces and coating, color change and visible residues. Damaged instruments must be disposed of (information on limitation of the numbers of processing cycles can be found in chapter 2.2). Instruments still having residues are to be cleaned and disinfected again (refer to chapter 5 above).

**8 Maintenance**

Instrument oils or grease must not be applied. Other maintenance aspects after application of chapters 5, 6 and 7 above are not necessary.

**9 Packaging for sterilization**

The cleaned and disinfected instruments have to be placed in the corresponding sterilization trays, which are then to be packed in single see-through sterilization pouches, that fulfil the following requirements (material/process):

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 142 °C (288 °F), sufficient steam permeability)

- sufficient protection of the instruments as well as of the sterilization packagings to mechanical damage

The packaging should be large enough to ensure that the seal is not under tension. In addition, do not use packaging that is too large to prevent the instruments from sliding around in the packaging.

A maximum weight of 0.26 kg per packaging is validated (not to be exceeded to remain within validated process).

**10 Sterilization**

Sterilization must be carried out using only the listed sterilization procedures.

Steam sterilization

- fractionated vacuum/dynamic air removal procedure<sup>4,5</sup> (with sufficient product drying<sup>6</sup>)
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665)
- sterilization time (exposure time at the sterilization temperature):

area	fractionated vacuum/ dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min <sup>6</sup>	not recommended
Germany	at least 5 min <sup>7</sup> at 134 °C (273 °F)	not recommended
other countries	at least 3 min <sup>7</sup> at 132 °C (270 °F) / 134 °C (273 °F)	not recommended

*The fundamental suitability of the instruments for an effective steam sterilization was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of a steam sterilizer compliant to EN ISO 17665 and EN ISO 13060) and the fractionated vacuum/dynamic air removal procedure. For this, typical conditions in clinic and doctor's practice as well as the specified procedure were considered.*

The flash/immediate use sterilization procedure must not be used.

Heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization must not be used.

<sup>4</sup> at least three vacuum steps

<sup>5</sup> The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure, requires significantly longer sterilization times as well as a sterilizer, procedure, parameter, and product specific validation under sole responsibility of the user.

<sup>6</sup> The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and by this is to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

<sup>7</sup> respectively 18 min (inactivation of prions, not relevant for USA)

Other sterilization procedures must not be used as processes are not validated and effects are not available.

### 11 Storage

The instruments, after sterilization, must be stored in the intact sterilization packaging at a dry and dust-free place.

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