

# Intensiv ApproxOpener

Serrated diamond coated metal strip for manual opening of tight contact points in Orthodontics

## Indications

- Manual approximal opening of tight contact points in orthodontic treatments (stripping)

## Benefits

- Frictionless application between contact points thanks to ultra-fine grain size and serrated edge
- Calibrated opening for application of the stripping method
- Safe application due to arch tension

Dr. Francesco Garino, Turin, Italy



## Product Description

Diamond coated metal strip clamped into a small metal bow; in 2 versions:

- Serrated and one-sided diamond coating (8 $\mu$ m), thickness approx. 0.09 mm
- Serrated and double-sided diamond coating (8 $\mu$ m), thickness approx. 0.10 mm
- Height of strip: 4 mm
- Length of strip: 20 mm



REF AO2018DS/3 Double Sided  
REF AO2018OS/3 One Sided