

Intensiv Composhape Set 1

Diamond instruments for modelling and finishing composite restorations

Contouring and finishing of composite fillings is an important step in restorative dentistry. Tooth substance must be preserved and enamel and filling fractures avoided to prevent infiltrations and secondary caries.



REF 015

Product Description

- 2 flames, Ø 010 and 012, 5 and 7 mm, 40 µm
- 1 torpedo, Ø 014, 6mm, 40 µm
- 1 football, Ø 021, 5mm, 40 µm
- 2 balls, Ø 017 and 024, 40 µm
- Shape and dimension matching instruments in 15 µm

Colour code	Micron	Description
514	40 µm	Fine
504	15 µm	Extra fine

LENGHT in mm L

COLOUR CODE Grit size indicator

DIAMETER in ISO 1/10 mm

REF 015 - FG	L 5.0 Ø 012	L 7.0 Ø 010	L 6.0 Ø 014	L 5.0 Ø 021	Ø 017	Ø 024
FG 4205	FG 4205L	FG 4062	FG 4255	FG 4201	FG 4400	All instruments are also available in RA
FG 5205	FG 5205L	FG 5062	FG 5255	FG 5201	FG 5400	

Indications

Contouring and finishing of:

- Composite fillings in anterior and posterior teeth
- Extended fissure sealings with flowable composites
- Enamel margins in adhesive preparation

Benefits

- Significant reduction of enamel and filling margin fractures
- Smooth filling surfaces ideally prepared for high gloss polishing

University of Zurich, Switzerland

