

Intensiv Geneva Prep Set

Diamond instruments for crown and bridge prosthetic restorations

Instruments for classic crown and bridge restorations according to methods developed at the University of Geneva.



REF 055

Product Description

- 1 cone, round, Ø 012, 7 mm, 90 µm
- 2 cones, round, Ø 016 and 018, 10 and 8 mm, 106 µm
- 2 cones, round, Ø 022 and 023, 10 and 11 mm, 106 µm
- 1 football Ø 023, 5 mm, 106 µm
- 2 cones, round, Ø 016 and 024, 6 and 7 mm, 50 µm
- 1 football, Ø 016, 3 mm, 50 µm

Colour code	Micron		Description
524	106 µm	60/80/90 µm*	Standard
514	50 µm		Golden Burs GB

*Depending on the shape and size of the instruments, the grit size may differ from the specified value.

Instrument	Length (L) in mm	Diameter (Ø) in mm	Material
FG D16	7.0	012	Standard
FG D6	10.0	016	Standard
FG 235	8.0	018	Standard
FG 237	10.0	022	Standard
FG 240	11.0	023	Standard
FG 255	5.0	023	Standard
FG D8 GB	6.0	016	Golden Burs GB
FG D7 GB	7.0	023	Golden Burs GB
FG 274 GB	3.0	016	Golden Burs GB

Indications

- Stump preparation for crowns and bridges

Benefits

- Reduced number of instruments for complex preparation procedures
- Method tested by the University of Geneva

University of Geneva, Switzerland

