

# Intensiv InstrumentTray

Compact, high quality stainless steel instrument tray for the storage and transport of long Intensiv instruments

During the processes in the dental practice, it's important to benefit from a support which is safe, hygienic, sterilizable and customizable in the content. The identification of the instrument to be used has to be clear and immediate. The instruments should be stored in absolute safety.



REF T600

## Product Description

- Tray in stainless steel, consisting in a base and lid
- The Tray has an internal folding support with six holes suitable to support sterilisable coloured silicone holders
- Dimensions: 87x51x14 mm (LxWxH)
- Maximum instrument length 35 mm
- The Tray has slots for the disinfection liquid to flow through it
- Coloured silicone holders allow the assignment of six Intensiv instruments

01 Coloured silicone holders in the tray allow allocation of the Intensiv instruments










02 Selection of Intensiv instruments



01



02

Bur-shank silicone holders									
REF	056green/10	056grey/10	056blu/10	056brown/10	056red/10	056white/10	056yellow/10	056orange/10	056ass/8
µm (the grit size is related to the strips)	125	90	80	60	40	25	15	08	

All coloured silicone holders are available in packaging of 10 pieces or in assortment of 8 pieces. 8 silicone holders are delivered free with the Tray. The desired silicone holders are to be ordered with the Tray.

## Application

- Safe storage and transport of Intensiv instruments in the dental practice

## Benefits

- The Intensiv instruments are available during the treatment according to the assignment of their use
- Coloured silicone supports for individual combination of the needed Intensiv instruments
- Compact and light
- Suitable for disinfection and sterilization
- Storage in complete safety of Intensiv instruments