

# Intensiv Mini Prep Set

Diamond instruments for minimally invasive restorations

Minimally invasive restoration techniques are designed to preserve healthy tooth substance. A defect-oriented procedure therefore requires preparation techniques with appropriately miniaturized instrument shapes.



REF 144

## Product Description

- 1 cylinder, Ø 007, 2.0 mm, 40 µm
- 2 balls, Ø 007, with various neck lengths, 40 µm
- 1 wheel, Ø 015, 0.7 mm, 40 µm

Colour code	Micron	Description
■ 514	40 µm	Fine

  
  

REF 144 - FG				
	■ FG 4612	■ FG 4199	■ FG 4699	■ FG 4610

## Indications

Preparation of:

- Class I amalgam and composite fillings
- Classes II and III composite fillings
- Palatal fillings
- Tunnel preparations

## Benefits

- Preservation of healthy tooth substance
- Preservation of marginal ridge in tunnel preparations
- Selective access
- No need of final polishing

University of Zurich, Switzerland



## Clinical pictures

- 01 Depth control of the preparation with instrument FG 4612
- 02 Preparation of the palatal micro cavity with instrument FG 4610
- 03 Palatal opening and preparation of an anterior cavity with instrument FG 4610