

Intensiv Ortho-StripsTray

Compact stainless-steel tray with coloured silicone supports for the placement of the Intensiv Ortho-Strips

During the processes in the dental practice, it's important to benefit from a support which is safe, hygienic, sterilizable and customizable in the content. The identification of the instrument to be used has to be clear and immediate. The storage of the strips should be carried out in complete safety.


















REF OST400

Product Description

- Tray in stainless steel, consisting in a base and lid
- The Tray has an internal folding support with six holes suitable to support sterilisable coloured silicone holders
- Dimensions: 87x51x14 mm (LxWxH)
- The Tray has slots for the disinfection liquid to flow through it
- Coloured silicone holders allow the arrangement of six different Intensiv Ortho-Strips

- 01 Introduction of the coloured silicone supports in the Tray for the positioning of the Intensiv Ortho-Strips
- 02 Selection of Intensiv Ortho-Strips



Bur-shank silicone holders							       
REF	056blu/10	056brown/10	056red/10	056white/10	056yellow/10	056orange/10	056ass/8
µm (the grit size is related to the strips)	80	60	40	25	15	08	

All coloured silicone holders are available in packaging of 10 pieces or in assortment of 8 pieces. 8 silicone holders are delivered free with the Tray. The desired silicone holders are to be ordered with the Tray.

Application

- Secure storage of Intensiv Ortho-Strips in the dental practice

Benefits

- Intensiv Ortho-Strips are assigned according to colour/grain
- Coloured silicone supports for individual combination
- Compact and light
- Suitable for disinfection and sterilization
- Storage and transport in complete safety of Intensiv Ortho-Strips