

Intensiv Perio Set

Diamond instruments for odontoplasty and mechanical root planing in periodontal treatments

Some important interventions in the treatment of periodontitis include the cleaning and polishing of the root surface. Supragingival and subgingival plaque as well as calculus and superficial endotoxin-containing cementum layers must be thoroughly removed. These are absolute preconditions for the complete healing and regeneration of the periodontal tissue.



REF 045, Intensiv Perio Set.
All Diamond instruments are available in packages of 1, 3 or 6 pcs

Product Description

- The coarse 75 µm burs are used only in odontoplasty
- The fine grit 40 µm burs are used for depuration of root surfaces
- The 15 µm burs are used for the final planing

Colour code	Micron	Description
524	80 µm	Medium
514	40 µm	Fine
504	15 µm	Extra fine

REF 045 – 831L – RA	REF 045 – 832L – RA
 RA 475	 RA 440
 RA 415	 RA 675
 RA 575	 RA 640
 RA 540	 RA 615
 RA 515	 RA 775
	 RA 740
	 RA 715

Indications

- Mechanical removal of supragingival and subgingival concretions
- Periodontal surgery (depuration of exposed root surfaces)
- Root planing
- Odontoplasty

Benefits

- Homogeneous smooth root surfaces

Prof. K.H. Rateitschak,
University of Basel, Switzerland

