

Intensiv ProxoshapeTray

Compact stainless steel tray with coloured silicone holders for the support of Intensiv Proxoshape oscillating files

During the processes in the dental practice, it is important to benefit from a support which is safe, hygienic, sterilizable and customizable in content. The identification of the instrument to be used has to be clear and immediate. The storage of the files should be carried out in complete safety.




REF PST500

Product Description

- Tray in stainless steel, consisting in a base and lid
- The Tray has an internal folding support with nine holes suitable to support sterilisable coloured silicone holders
- Dimensions: 87x51x14 mm (LxWxH)
- The Tray has slots for the disinfection liquid to flow through it
- Coloured silicone holders allow the arrangement of different oscillating Intensiv files

- 01 Easy identification of Intensiv files
- 02 Secure Intensiv file removal



Bur-shank silicone holders									
REF	056green/10	056grey/10	056blu/10	056brown/10	056red/10	056white/10	056yellow/10	056orange/10	056ass/8
µm (the grit size is related to the strips)	125	90	80	60	40	25	15	08	

All coloured silicone holders are available in packaging of 10 pieces or in assortment of 8 pieces. 9 silicone holders as shown in the picture are delivered free with the Tray. The desired silicone holders are to be ordered with the Tray

Application

- Secure storage of Intensiv oscillating files in the dental practice

Benefits

- Intensiv files are assigned according to colour/grain
- Coloured silicone supports for individual combination
- Compact and light
- Suitable for disinfection and sterilization
- Storage and transport in complete safety of Intensiv files