

Intensiv Sets Direct and Indirect Anterior & Posterior

Diamond instruments for anterior and posterior direct and indirect restorations

Aesthetic direct and indirect restorations play a predominant role in modern dentistry and represent a new limit between the conservative dentistry and the prosthesis for recovering individual elements. A complete set of instruments for all clinical indications in the field of direct and indirect aesthetic restorations, is selected.



REF 202PI



REF 202PD



REF 202ADI

Product Description

3 sets with 8 diamond-coated instruments each, placed in the Intensiv HygienicTray HT1000.

REF 202PI

Posterior indirect restorations

- 1 cone, Ø 012, 8 mm, 80 µm
- 4 cones, Ø 014, 8 mm, 80 and 25 µm grits
- 1 cone, Ø 020, 4 mm, 25µm
- 2 cones, Ø 018, 6 mm, 80 and 25 µmm

REF 202PD

Posterior direct restorations

- 1 cylinder, Ø 009, 3 mm, 80 µm
- 4 cylinders, Ø 011 and 013, 6 mm, 80 and 25 µm
- 2 cylinders, Ø 014, 6 mm, 80 and 25 µm

REF 202ADI

Anterior direct and indirect restorations

- 1 ball, Ø 009, 60 µm
- 2 balls, Ø 014, 90 and 40 µm
- 2 balls, Ø 023, 106 and 40 µm
- 2 cylinders, rounded end, Ø 012, 10 mm, 90 and 40 µm
- 1 flame, Ø 010, 7 mm, 40 µm

Indications

- Preparation of cavities for direct and indirect restorations in anterior and posterior areas
- Conservative preparations with perfectly finished cavity walls and margins

Benefits





- Placed in Intensiv HygienicTray (HT1000) instrument holders, to facilitate the identification and the organization of the sets
- Possible individualisation of the sets with additional instrument insertion (oscillating instruments, rubber cups etc.) thanks to the 13 FG and 5 RA holes in the tray
- Form and size congruence between the preparation and finishing instruments

Dr. Roberto Spreafico, Busto Arsizio, Italy




















Clinical pictures

- 01 FG 8525: preparation of the walls
- 02 FG 3413R: finishing of the occlusal surface
- 03 FG 8212S: opening of the cavity for direct restorations
- 04 FG 3514/3614B/3414B: finishing of the margins
- 05 FG 305L: production of vestibular guiding grooves and preparation of incisal margin
- 06 FG 4307N: finishing of the preparation for veneers

Colour code	Micron		Description
 524	106 µm	60/80/90 µm*	Standard
 524	80 µm		Medium
 514	40 µm		Fine
 514	25 µm		Fine

*Depending on the shape and size of the instruments, the grit size may differ from the specified value.



REF 202PI – FG								
REF 202PD – FG								
REF 202ADI – FG	