

Intensiv Tooth Caring Bur Set

Diamond instruments for lowering and finishing subgingival edge of cervical steps and shoulders of fixed prosthetic preparations without damaging adjacent teeth

This preparation method requires a set of special instruments to enable efficient and controlled substance removal. At the same time these instruments ensure perfectly finished surfaces of shoulders and margins for stumps in prosthetics.



REF 001

Product Description

- 2 cylinders, rounded edges, end-coated only, Ø 011, 40 and 25 µm
- 3 shape and dimension matching instruments, Ø 014 in 80, 40 and 25 µm

Colour code	Micron	Description
524	80 µm	Medium
514	40 µm	Fine
514	25 µm	Fine

REF 001 – VB1 – FG	Ø 011	Ø 011	REF 001 – VB2 – FG	Ø 014	Ø 014	Ø 014
	FG 01140	FG 01125		FG 01480	FG 01440	FG 01425

Indications

- Preparation of shoulder margins and proximal box margins

Benefits

- Preservation of prepared stump and shoulder profiles
- Prevention of iatrogenic damage

University of Bern, Switzerland

