

Intensiv Universal full Crown and Porcelain Veneer Prep Set

Diamond Instruments for crown and bridge prosthetic restorations

Often during prosthetic preparation it is necessary to follow a simple, accurate, repeatable and reliable method. The method created by Dr. Dazhaev, Moscow, in collaboration with Intensiv, allows for a structured process, thanks to the adequate selection of instruments and the production of special shapes, ensuring long-term outcome.



Product Description

- 3 cylinders, round, Ø 010, 012 and 021, 8 mm, 125 µm
- 2 cylinders, round, Ø 013 and 022, 8 and 5 mm, 40 µm
- 2 egg-shaped instruments, Ø 023 and 022, 5 mm, 125 and 40 µm
- 2 cones, round, Ø 012 and 013, 10 mm, 125 and 40 µm
- 4 depth markers, 125 µm

REF 010

Colour code	Micron		Description
534	125 µm	106 µm*	Coarse
514	40 µm		Fine

*Depending on the shape and size of the instruments, the grit size may differ from the specified value.

Instrument	Length (L) in mm	Diameter (Ø) in mm	FG Code
REF 010 - FG	8.0	010	FG 309 C
	8.0	012	FG 321 C
	8.0	013	FG 4315S
	8.0	021	FG 235A C
	8.0	022	FG 4235S
	5.0	023	FG 250 C
	10.0	012	FG 305L C
	10.0	013	FG 4307N
REF 010 - FG	5.0	022	FG 4250B
Depth - Marker	2.0	010	FG 707 C
	1.0	Penetration depth 004	FG 102A C
	1.0	Penetration depth 006	FG 103A C
	1.0	Penetration depth 004	FG 108A C

Indications

- Preparation of full crowns
- Preparation of porcelain veneers

Benefits

- Fast and very precise reduction
- Depth markers
- Very accurate crisp chamfer creation
- Composition of this set allows following both, classic and minimal invasive approach to preparation

Dr. G. Dazhaev, Moscow, Russia

